

PURE
COMFORT
PURE
PERFORMANCE

WHEELCHAIR
CUSHIONS



Elasto-GelTM



southwest technologies inc.

"Treating the world well"®



southwest technologies inc.

1746 Levee Road, North Kansas City, MO 64116
ph:(800) 247-9951 ph:(816) 221-2442 fax:(816) 221-3995
e-mail: info@elastogel.com • website: www.elastogel.com

Distributed by:

To:

Quality Policy

It is the policy of Southwest Technologies Inc. to design, manufacture, market and distribute products consistent with our high quality standards. We periodically review the products for design and performance improvements.


The products will be designed and manufactured to reliably perform their intended functions as labeled. All products will be produced and marketed within the established company guidelines.

Every component of the



is designed to assure comfort & performance.



- Unique triple foam design assures consistent body positioning.
- Soft Elasto-Gel™ combined with special resilient foams allows efficient weight disbursement off ischials and coccyx.
- Elasto-Gel™ component helps maintain a cool temperature and has vapor transmission properties to help keep the body dry and comfortable.
- Available in sizes to fit most wheelchair needs and body weights.
- Optional insert is available to overcome the negative postural effects of sling seats.






Warranty Information: Southwest Technologies' (SWT) wheelchair cushions The Elasto-Gel™ cushions are warranted to be free of defective workmanship and materials at the time of purchase and are guaranteed for eighteen (18) months.

Positioning Contour Cushion

- High Glycerine Gel-High Memory
- Positioning Foam
- Top Foam
- Bottom Foam



SC4021 Firmer Foam
Recommended for patients weighing 100lbs. to 185lbs.



SC6021 Soft & Firmer Foam
Recommended for patients weighing 185lbs. to 280lbs.

Ordering Information

Stock #	Description	UPC	Packaging
SC4021 (A)	16"x18" Positioning/Contour	745713-04020	1 each
SC4021 (B)	17"x17" Positioning/Contour	745713-04021	1 each
SC4021 (C)	18"x16" Positioning/Contour	745713-04022	1 each
SC6021 (A)	16"x18" Positioning/Contour	745713-05023	1 each
SC6021 (B)	17"x17" Positioning/Contour	745713-05024	1 each
SC6021 (C)	18"x16" Positioning/Contour	745713-05025	1 each
SC2100	Leveling Insert	745713-02100	1 each

Always specify size when ordering.

HCPC Code: E2603